EL sensor newsletter

NÚMERO 1:2021

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EL SENSOR: Newsletter, Número 1

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1. Dover AR et al. Assessment of the effect of the COVID-19 lockdown on glycaemic control in people with type 1 diabetes using flash glucose monitoring. *Diabet Med*. 2020 Aug 2. doi: 10.1111/dme.14374
2. Fernandez E et al. Impact of COVID-19 lockdown on glycaemic control in patients with type 1 diabetes. *Diabetes Res Clin Pract* 2020;166 doi: 10.1016/j.diabres.2020.108348
3. Prabhu Navis J et al. Impact of COVID-19 lockdown on flash and real-time glucose sensor users with type 1 diabetes in England. *Acta Diabetol*. 2020 Oct 16:1-7. doi: 10.1007/s00592-020-01614-5
4. Boscardi F, et al. Glycaemic control among people with type 1 diabetes during lockdown against the SARS-CoV-2 outbreak in Italy. EASD 2020, Oral Presentation # 104

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1. Alva S, et al. Accuracy of a 14-Day Factory-Calibrated Continuous Glucose Monitoring System With Advanced Algorithm in Pediatric and Adult Population With Diabetes. *J Diabetes Sci Technol*. 2020 Sep 19;1932296820958754. doi: 10.1177/1932296820958754
2. Bailey T, et al. The performance and usability of a factory-calibrated flash glucose monitoring system. *Diabetes Technol Ther*. 2015;17(11):787-794.

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Wilmot EG et al. Time in range: A best practice guide for UK diabetes healthcare professionals in the context of the COVID-19 global pandemic. *Diabet Med*. 2020 Oct 19:e14433. doi: 10.1111/dme.14433

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Lu J, et al. Time in Range Is Associated with Carotid Intima-Media Thickness in Type 2 Diabetes Technol Ther. 2020 Feb;22(2):72-78. doi: 10.1089/dia.2019.0251. Epub 2019 Oct 11

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Seibold A. Minimising adverse skin reactions to wearable continuous glucose monitoring sensors in patients with diabetes. J Diabetes Sci Technol. 2021; Jan 8, doi: 10.1177/1932296820984763

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Okada J, Osaki A, Shimoda Y, et al. Differences in the effects of Kenyan, Tanzanian, and Ethiopian coffee intake on interstitial glucose levels measured by FreeStyle Libre: A pilot case study. Curr Ther Res Clin Exp. 2020 Sep 19;93:100606. doi: 10.1016/j.curtheres.2020.100606

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Yang J, et al. Association of time in range, as assessed by continuous glucose monitoring, with painful diabetic polyneuropathy. J Diabetes Investig. 2020 Sep 3. doi: 10.1111/jdi.13394

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Suzuki J et al. Association between scanning frequency of flash glucose monitoring and continuous glucose monitoring-derived glycemic makers in children and adolescents with type 1 diabetes. Pediatr Int. 2020 Aug 3. doi: 10.1111/ped.14412

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Al Hayek AA, et al. Effectiveness of the FreeStyle Libre Flash Glucose Monitoring System on Diabetes Distress Among Individuals with Type 1 Diabetes: A Prospective Study. *Diabetes Ther*. 2020 Apr;11(4):927-937

Al Hayek AA, et al. Assessing Diabetes Distress and Sleep Quality in Young Adults with Type 1 Diabetes Using FreeStyle Libre: A Prospective Cohort Study. *Diabetes Ther*. 2020 Jul;11(7):1551-1562. doi: 10.1007/s13300-020-00849-3

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Stimson RH, et al. HbA1c Is Disproportionately Higher in Women and Older People With Type 1 Diabetes Compared With Flash Glucose Monitoring Metrics of Glycemic Control *J Diabetes Sci Technol*. 2020 Oct 26; doi:10.1177/1932296820967335

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Pease A et al. Cost-effectiveness of health technologies in adults with type 1 diabetes: a systematic review and narrative synthesis. *Systematic Reviews* 2020; 9: 171 <https://doi.org/10.1186/s13643-020-01373-y>

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Biester T et al. Intermittently scanned glucose values for continuous monitoring: Cross-sectional analysis of glycaemic control and hypoglycaemia in 1809 children and adolescents with type 1 diabetes. *Diabetes Technol Ther*. 2020 Oct 21. doi: 10.1089/dia.2020.0373

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Moser O, et al. Glucose management for exercise using continuous glucose monitoring (CGM) and intermittently scanned CGM (isCGM) systems in type 1 diabetes: position statement of the European Association for the Study of Diabetes (EASD) and of the International Society for Pediatric and Adolescent Diabetes (ISPAD) endorsed by JDRF and supported by the American Diabetes Association (ADA). *Pediatr Diabetes*. 2020; 21(8):1375-1393. doi: 10.1111/peci.13105.

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Gal RL, et al. Diabetes Telehealth Solutions: Improving Self-Management Through Remote Initiation of Continuous Glucose Monitoring. *J Endocr Soc*. 2020 Jun 23;4(9):bvaa076. doi: 10.1210/endo/bvaa076

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1. Battelino T, Danne T, Bergenstal RM, Amiel SA, Beck R, Biester T, et al. Clinical Targets for Continuous Glucose Monitoring Data Interpretation: Recommendations From the International Consensus on Time in Range. *Diabetes Care* 2019; 42(8): 1593-1603 doi:10.2337/dci19-0028

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Department of Pediatric Endocrinology, Diabetes and Metabolism, University Children's Hospital, University Medical Center (Ulm), and Faculty of Medicine, University of Ulm, Germany
 2Diabetes Centre for Children and Adolescents, Kinder- und Jugendklinikhaus Auf der Suut, Hannover, Germany
 3International Diabetes Center at Park College, Milledgeville, MN
 4Diabetes Research Group, Kings College London, London, UK
 5Abbott Center for Health Research, Tampa, FL
 6Children Research Institute, MC25 Sun Health Hospital, Vita Salute San Raffaele University, Milan, Italy
 7Division of Endocrinology and Diabetes, Department of Pediatrics, Stanford Medical Center, Stanford, CA

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1Department of Pediatric Endocrinology, Diabetes and Metabolism, University Children's Hospital, University Medical Center (Ulm), and Faculty of Medicine, University of Ulm, Germany
 2Diabetes Centre for Children and Adolescents, Kinder- und Jugendklinikhaus Auf der Suut, Hannover, Germany
 3International Diabetes Center at Park College, Milledgeville, MN
 4Diabetes Research Group, Kings College London, London, UK
 5Abbott Center for Health Research, Tampa, FL
 6Children Research Institute, MC25 Sun Health Hospital, Vita Salute San Raffaele University, Milan, Italy
 7Division of Endocrinology and Diabetes, Department of Pediatrics, Stanford Medical Center, Stanford, CA
 8American Diabetes Association, Alexandria, VA

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1. Cho J-H, et al. Long-Term Effect of the Internet-Based Glucose Monitoring System on HbA1c Reduction and Glucose Stability: A 30-month follow-up study for diabetes management with a ubiquitous medical care system. *Diabetes Care*. 2006;29(12):2625-2631
2. Verhoeven F, et al. The Contribution of Teleconsultation and Videoconferencing to Diabetes Care: A Systematic Literature Review. *J Med Internet Res*. 2007;9(5):e37
3. Charpentier G, et al. The Diabeo Software Enabling Individualized Insulin Dose Adjustments Combined With Telemedicine Support Improves HbA1c in Poorly Controlled Type 1 Diabetic Patients. *Diabetes Care*. 2011;34(3):533-539
4. Dixon RF, et al. A Virtual Type 2 Diabetes Clinic Using Continuous Glucose Monitoring and Endocrinology Visits. *J Diabetes Sci Technol*. 2020;14(5):908-911
5. Bergenstal RM, et al. Remote Application and Use of Real-Time Continuous Glucose Monitoring by Adults with Type 2 Diabetes in a Virtual Diabetes Clinic. *Diabetes Technol Ther* 2020; doi: 10.1089/dia.2020.0396
6. Polonsky WH, et al. Impact of Participation in a Virtual Diabetes Clinic on Diabetes-Related Distress in Individuals With Type 2 Diabetes. *Clin Diabetes*. 2020;38(4):357-362
7. Wood CL, et al. Use of Telemedicine to Improve Adherence to American Diabetes Association Standards in Pediatric Type 1 Diabetes. *Diabetes Technol Ther*. 2016;18(1):7-14
8. von Sengbusch, et al. Outcomes of monthly video consultations as an add-on to regular care for children with type 1 diabetes: A 6-month quasi-randomized clinical trial followed by an extension phase. *Pediatr Diabetes*. 2020;21(8):1502-1515
9. von Sengbusch S, et al. Parental expectations before and after 12-month experience with video consultations combined with regular outpatient care for children with type 1 diabetes: a qualitative study. *Diabetic Med*. 2020;e14410
10. Gal RL, et al. Diabetes Telehealth Solutions: Improving Self-Management through Remote Initiation of Continuous Glucose Monitoring. *J Endocr Soc*. 2020;4(9):bvaa076. doi: 10.1210/endo/bvaa076
11. Vigersky RA, et al. The Effectiveness of Virtual Training on the MiniMed 670G System in People with Type 1 Diabetes During the COVID-19 Pandemic. *Diabetes Technol Ther* 2020; doi: 10.1089/dia.2020.0234

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1. Hermanns N, et al. The impact of a structured education and treatment programme (FLASH) for people with diabetes using a flash sensor-based glucose monitoring system: Results of a randomized controlled trial. *Diabetes Res Clin Pract.* 2019 ; 150:111-121
2. Schluter S, et al. Evaluation of the SPECTRUM training programme for real-time continuous glucose monitoring: A real-world multicentre prospective study in 120 adults with type 1 diabetes. *Diabet Med.* 2020 Nov 23;e14467. doi: 10.1111/dme.14467. Online ahead of print.
3. Pemberton JS, et al. DYNAMIC: DYNAMIC glucose Management strategies delivered through a structured education program improves time in range in a socioeconomically deprived cohort of Children and young people with type 1 diabetes with a history of hypoglycaemia. *Pediatr Diabetes* 2020; Nov: doi: 10.1111/pedi.13155

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1. National Diabetes Statistics Report 2020, Centers for Disease Control and Prevention (CDC)
2. Diabetes Canada. https://diabetes.ca/DiabetesCanadaWebsite/media/Advocacy-and-Policy/Backgrounder/2020_Backgrounder_Canada_English_FINAL.pdf