




FreeStyle
Libre 2

Case study: Tony



Case study

Images are for illustrative purposes only. Not actual patient or data. The FreeStyle LibreLink app is only compatible with certain mobile devices and operating systems. Please check the website for more information about device compatibility before using the app. Use of FreeStyle LibreLink requires registration with LibreView. © 2023 Abbott. FreeStyle, Libre, and related brand marks are owned by Abbott. ADC-70040 v1.0 03/23.

 **Abbott**
life. to the fullest.®



Case study: Tony

Age	53	BMI	28.4 kg/m ²
Diabetes (Type)	Type 2 diabetes	Last HbA1c value	56 mmol/mol (7.3%)
Profession	Machine operator (night shift)	Target glucose range	3.9–10 mmol/L
Duration of diabetes	6 years	Treatment	Basal-bolus insulin therapy



Summary

Tony is a machine operator at a local factory, who developed Type 2 diabetes 6 years ago. His job is quite physical, and he regularly works the night shift.



Specific objective

Reduce the incidence of hypoglycaemia late in the evening and overnight while Tony is working.

Case study: Tony



AGP Report

11 May 2021 - 24 May 2021 (14 Days)


GLUCOSE STATISTICS AND TARGETS

11 May 2021 - 24 May 2021 **14 Days**
 % Time Sensor is Active **89%**

Ranges And Targets For		Type 1 or Type 2 Diabetes
Glucose Ranges	Targets	% of Readings (Time/Day)
Target Range 3.9-10.0 mmol/L	Greater than 70% (16h 48min)	
Below 3.9 mmol/L	Less than 4% (58min)	
Below 3.0 mmol/L	Less than 1% (14min)	
Above 10.0 mmol/L	Less than 25% (6h)	
Above 13.9 mmol/L	Less than 5% (1h 12min)	
Each 5% increase in time in range (3.9-10.0 mmol/L) is clinically beneficial.		


Average Glucose **7.4 mmol/L**
Glucose Management Indicator (GMI) **6.3% or 45 mmol/mol**
Glucose Variability **32.3%**
 Defined as percent coefficient of variation (%CV); target ≤36%

AMBULATORY GLUCOSE PROFILE (AGP)



LibreView

TIME IN RANGES



Snapshot


11 May 2021 - 24 May 2021 (14 Days)

LibreView

Glucose


GMI **6.3 %** or **45 mmol/mol**

AVERAGE GLUCOSE **7.4 mmol/L**
 % above target **23 %**
 % in target **69 %**
 % below target **8 %**




LOW GLUCOSE EVENTS

15
 Average duration **110 Min**



Sensor Usage

% TIME SENSOR IS ACTIVE **89 %**
 Daily scans **9 / Day**




Images are for illustrative purposes only. Not actual patient data.

Case study: Tony


Weekly Summary LibreView

11 May 2021 - 24 May 2021 (14 Days)



Weekly Summary LibreView

11 May 2021 - 24 May 2021 (14 Days)



Images are for illustrative purposes only. Not actual patient data.

What does the 4-step review tell us?

STEP 1

Data capture and Time in Range (TIR)

Tony's 89% data capture and his Time in Range of 69% are very good achievements. He should be commended for this, especially with his work pattern.

STEP 2

Look for patterns of hypoglycaemia

There is a significant risk of hypoglycaemia in the evenings and overnight; Tony's median line skirts the lower limit of his target range in the evening, with his blue and grey bands of variability straying below this level from 8:00pm. A look at Tony's Snapshot report shows that he is experiencing low glucose events in the later evening and overnight, often below 3.0 mmol/L. Given that this is when Tony is at work, this requires urgent action.

STEP 3

Look for patterns of hyperglycaemia

There is a gradual rise in the blue median line from about 2:00am through breakfast time. It then shows another peak from 3:00pm until 7:00pm. These elevations are matched by Tony's blue shaded band, suggesting a consistent trend, as confirmed in the daily profiles in his Weekly Summary reports. However, his Time in Range is currently 69% and he spends 23% of time above the 10 mmol/L threshold, so there is nothing to be concerned about here.

STEP 4

Look for patterns of glucose variability

The width of Tony's blue and grey bands indicate low variability, except overnight where the grey band is wider. Again, this could be reflective of Tony's shift patterns and his CV is below 36%, indicating a stable glucose variability profile.

What actions might you agree with Tony?

- A reduction in Tony's basal insulin is recommended to reduce the evening and overnight incidence and risk of hypoglycaemia.
- It would be beneficial to understand more about the demands of Tony's job and shift patterns to further tailor his medication.

Case study: Tony

AGP Report

4 August 2021 - 17 August 2021 (14 Days)

GLUCOSE STATISTICS AND TARGETS

4 August 2021 - 17 August 2021 **14 Days**
% Time Sensor is Active **93%**

Ranges And Targets For Type 1 or Type 2 Diabetes

Glucose Ranges	Targets % of Readings (Time/Day)
Target Range 3.9-10.0 mmol/L	Greater than 70% (16h 48min)
Below 3.9 mmol/L	Less than 4% (58min)
Below 3.0 mmol/L	Less than 1% (14min)
Above 10.0 mmol/L	Less than 25% (6h)
Above 13.9 mmol/L	Less than 5% (1h 12min)
Each 5% increase in time in range (3.9-10.0 mmol/L) is clinically beneficial.	


Average Glucose **7.8** mmol/L

Glucose Management Indicator (GMI) **6.5% or 48** mmol/mol

Glucose Variability **33.8%** mol


Defined as percent coefficient of variation (%CV); target ≤36%

AMBULATORY GLUCOSE PROFILE (AGP)



LibreView

TIME IN RANGES



Snapshot


4 August 2021 - 17 August 2021 (14 Days)

LibreView


Glucose

GMI **6.5%** or **48** mmol/mol

AVERAGE GLUCOSE **7.8** mmol/L
 % above target **33%**
 % in target **62%**
 % below target **4%**




LOW GLUCOSE EVENTS **5**
 Average duration **132** min



Sensor Usage

% TIME SENSOR IS ACTIVE **93%**
 Daily scans **5 / Day**




Images are for illustrative purposes only. Not actual patient data.

Case study: Tony


Weekly Summary *LibreView*

4 August 2021 - 17 August 2021 (14 Days)



Weekly Summary *LibreView*

4 August 2021 - 17 August 2021 (14 Days)



Images are for illustrative purposes only. Not actual patient data.

What does the 4-step review tell us?

STEP 1

Data capture and Time in Range (TIR)

Tony's data capture is excellent at 93%, however his Time in Range has declined to 62% from 69%, so it is worth understanding why.

STEP 2

Look for patterns of hypoglycaemia

There has been an improvement in Tony's Time Below Range, both in the amount of time below 3.9 mmol/L and below 3.0 mmol/L. Tony's **Snapshot report** shows that he is experiencing fewer low glucose events, but that he is at risk of hypoglycaemia late in the evening and overnight, from 9:00pm onwards.

STEP 3

Look for patterns of hyperglycaemia

There is still undulation in Tony's AGP profile around mealtimes, that takes him above 10 mmol/L, such that his Time Above Range has increased from 23% to 33%, with a reduction in Time in Range from 69% to 62%. This is confirmed again the daily glucose profiles in his **Weekly Summary report** and makes this a focus for improvement.

STEP 4

Look for patterns of glucose variability

Tony's blue and grey bands continue to indicate low variability, especially during the day. His grey band is wider from 8:00pm and through the night. Although his CV has increased to 33.8%, it is still below 36%.

What actions might you agree with Tony?

- Tony should be encouraged to lose weight; with his duration of type 2 diabetes this may provide an opportunity for the elimination of insulin to further reduce the risk of hypoglycaemia.
- Given his current lifestyle with varying and unpredictable physical activities, Tony should be congratulated for maintaining his current glucose control.