



FreeStyle Libre

FreeDM2 Randomised Controlled Trial (RCT) Results.

CONTENTS:

Global challenge	2
A first-of-its-kind RCT in T2D	3
Study design	4
HbA1c at 16 and 32 weeks	5
TIR at 16 and 32 weeks	6
User-driven changes	7
Behavioural change and long-term benefit	8

RCT=randomised controlled trial; T2D=type 2 diabetes; TIR=Time in Range.
Images are for illustrative purposes only. Not actual patient.



Abbott

Most basal insulin-treated T2D remains uncontrolled,¹ pointing to the need for more actionable data.

In the UK,
~70%¹
of PWD treated
with basal insulin
have an HbA1c
≥58 mmol/mol




Use of Libre systems is associated with an average HbA1c reduction of 0.32–1.4% in populations treated with basal insulin^{2–5}

PWD=people with diabetes; T2D=type 2 diabetes.

1. Interface audit data analysis, data held on file at Abbott. 2. Elliott, T. *Diabetes Vasc Dis Res* (2021): <https://doi.org/10.1177/14791641211021374>. 3. Carlson, A. L. *BMJ Open Diabetes Res Care* (2022): <https://doi.org/10.1136/bmjdr-2021-002590>. 4. Wright, E. E. *Diabetes Spectr* (2021): <https://doi.org/10.2337/ds20-0069>. 5. Nathanson, D. *Diabetologia* (2025): <https://doi.org/10.1007/s00125-024-06289-z>.

FreeDM2: unboxing the behavioural impact of using the Libre 3 system.

A rigorous landmark RCT to evaluate glycaemic and behavioural improvements in adults with T2D on basal insulin regimens with suboptimal HbA1c and TIR¹



A key strength of this study lies in its objective of understanding the mechanisms by which the Libre 3 system influences change in glycaemia without therapy intervention¹

RCT=randomised controlled trial; T2D=type 2 diabetes; TIR=Time in Range.

Images are for illustration only. Not actual patients or data.

1. Wilmot, E. G. *Lancet Diabetes Endocrinol* (2026): [https://doi.org/10.1016/S2213-8587\(26\)00076-8](https://doi.org/10.1016/S2213-8587(26)00076-8).

FreeDM2 RCT study design and endpoints.¹




Methodology:

- 24 centres
- 2:1 randomisation

Population:

- Mean age: 60.7 ± 9.8
- HbA1c: 7.5–11.0% (59–97 mmol/mol)
- Mean HbA1c: 8.8% (73 mmol/mol)
- Mean daily finger pricks: ~2
- T2D on basal insulin plus SGLT2i or GLP-1 or dual GIP/GLP-1 RA



PRIMARY ENDPOINT

16 HbA1c change at 16 weeks

SECONDARY ENDPOINT

32 HbA1c change at 32 weeks

CGM CGM metrics

Accelerometer Accelerometer data

Insulin Changes to insulin doses & non-insulin medications

Interviews Qualitative interviews

Safety Safety

PROMs PROMs (including EQ-5D-5L, GMSS, UKDDQ, HCS, HFS-II)

CGM=continuous glucose monitoring; EQ-5D-5L=EuroQoL 5-Dimension 5-Level questionnaire; GIP=glucose-dependent insulinotropic polypeptide; GLP-1 RA=glucagon-like peptide-1 receptor agonist; GMSS=Glucose Monitoring Satisfaction Survey; HCS=Hypoglycaemic Confidence Scale; HFS-II=Hypoglycaemia Fear Survey-II; PROMs=patient-reported outcome measures; SGLT2i=sodium-glucose cotransporter-2 inhibitor; SMBG=self-monitoring of blood glucose; UKDDQ=UK Diabetes and Diet Questionnaire.

1. Wilmot, E. G. *Lancet Diabetes Endocrinol* (2026); [https://doi.org/10.1016/S2213-8587\(26\)00076-8](https://doi.org/10.1016/S2213-8587(26)00076-8).

The use of the Libre 3 system resulted in clinically and statistically significant HbA1c changes at 16 weeks, sustained through 32 weeks.¹



Patients achieving $\geq 0.5\%$ reduction in HbA1c

16 weeks:

Libre 3 system	SMBG
65%	27%

$p < 0.0001$

32 weeks:

Libre 3 system	SMBG
70%	49%

$p = 0.0006$


SMBG=self-monitoring of blood glucose.

1. Wilmot, E. G. *Lancet Diabetes Endocrinol* (2026): [https://doi.org/10.1016/S2213-8587\(26\)00076-8](https://doi.org/10.1016/S2213-8587(26)00076-8). 2. Baseline HbA1c for study population was 8.81%.

The Libre 3 system users experienced increased Time in Range (TIR).¹



TIR improvements occurred with no increase in hypoglycaemia¹



T2D=type 2 diabetes.

1. Wilmot, E. G. *Lancet Diabetes Endocrinol* (2026): [https://doi.org/10.1016/S2213-8587\(26\)00076-8](https://doi.org/10.1016/S2213-8587(26)00076-8).

Real-time biofeedback helps people with T2D make behaviour changes that drive long-term glycaemic benefits.¹



T2D=type 2 diabetes.

Images are for illustration only. Not actual patients or data.

1. Wilmot, E. G. *Lancet Diabetes Endocrinol* (2026): [https://doi.org/10.1016/S2213-8587\(26\)00076-8](https://doi.org/10.1016/S2213-8587(26)00076-8). 2. This effect was recorded at 16 weeks (Phase A, self-management); no statistically significant difference at 32 weeks (Phase B, clinician-supported care).

FreeDM2 demonstrated improvement in glycaemic management and behavioural changes in adults with T2D on basal insulin regimens.¹



Most glycaemic improvements were achieved by patient behaviour changes rather than by medication changes and provider interventions¹



Significant reduction in HbA1c

- 1.0% for Libre 3 system at 32 weeks vs. 0.5% for SMBG at 32 weeks



Significant increase in TIR

- +2.5 hour total TIR improvement at 16 and 32 weeks



Proactive behavioural changes

- Healthier dietary choices²
- Increased light physical activity²



Psychosocial improvements

- Improved hypoglycaemia confidence
- Reduced emotional burden

SMBG=self-monitoring of blood glucose; T2D=type 2 diabetes; TIR=Time in Range.

Target glucose range: 3.9–10.0 mmol/L.

1. Wilmot, E. G. *Lancet Diabetes Endocrinol* (2026): [https://doi.org/10.1016/S2213-8587\(26\)00076-8](https://doi.org/10.1016/S2213-8587(26)00076-8). **2.** This effect was recorded at 16 weeks (Phase A, self-management); no statistically significant difference at 32 weeks (Phase B, clinician-supported care).

© 2026 Abbott. All Rights Reserved. Libre, the butterfly logo, the sensor shape and appearance, the colour yellow, and related marks and/or designs are the intellectual property of the Abbott group of companies in various territories. Other marks are the property of their respective owners. ADC-2689696 v2.0 05/26.

